

Cable guide

Fitting with Appraise 2.4	 17292002 (Airlink 2, GN Hearing) 17292000 (Airlink, GN Hearing)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 3122 150 73781 (FLEX CONNECTOR CLIP CS63) 19181700 (STRIP,FLEX,CS63,PRG)	 15206700 (Programming cable CS63, Right) 15206701 (Programming cable CS63, Left) 19181700 (STRIP,FLEX,CS63,PRG)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 19533900 (ADAPTER FLEX/3 PIN/CS53)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 18041400 (PROGRAMMING ADAPTOR,8PIN)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 17598200 (PRGM SHOE ASM,W/ CS44,DELHI)
Interton Centro						
RIE WL - CTx62	•	•	•	•		
Interton Cosmo						
RIE - Cx62		•	•	•		
RIE WL - Cx61	•	•	•	•		
Custom WL	•	•	•	•		
Interton Step						
BTE WL - STx80, STx70, STx65	•				•	
RIE WL - STx60	•				•	
Custom WL	•	•	•	•		
Custom		•	•	•		
Interton Crisp						
BTE WL - Clx80, Clx65	•				•	
RIE - Clx62		•	•	•		
RIE WL - Clx60	•				•	
Custom WL	•	•	•	•		
Custom		•	•	•		
Interton Gain						
BTE - GANx90						•
Interton Share						
BTE - SR1x80, SR1x70, SR1x65						•
RIE - SR1x60						•
Customs		•	•	•		

Please note that all Interton hearing instruments wireless and non-wireless **fittings require a fresh battery when fitted with Appraise 2.4.**